






Fine Particle Filter Installation Instructions

en

Instructions For Use

1 ➔ Remove the filter from the package.	
2 ➔ Place the filter into the housing with the foam (dark) side facing toward the ventilator. ➔ Ensure that the filter is flush with the inside surface of the housing.	
3 ➔ Installation complete.	

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.



Single use



Caution, consult accompanying documents

Part No. 10072519 Rev A

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.

© 2011 Covidien.
 Covidien Inc, 15 Hampshire Street,
 Mansfield, MA 02048 USA

www.covidien.com